

**Replacing the Bulb**

To replace the bulb (use only a 60W bulb), first unplug the clock from the wall. Remove the globe from the clock by pulling upward. Remove the globe by holding the base with one hand, and pulling up on the globe. Take care not to break the bulb by hitting it with the globe as you do this. **Caution: If the bulb was recently in use, it will be hot - let the lamp cool before removing the bulb.** Screw in the new bulb completely. Carefully put the globe back on by lowering it over the bulb (the slots on the globe face backwards) and make sure it is fully home and securely held. Plug in the clock, re-set the time, alarm time, and any other settings. Replacement bulbs are available from Bio-Brite or local hardware stores (ask for a 110 Volt 60 Watt E12 Candelabra base bulb).

**Warranty**

EZ Wake SunRise Clock is guaranteed against any defect in its manufacture for 12 months from the date of purchase.

**Service**

If your EZ Wake SunRise Clock requires service, please send it to: Bio-Brite Service Department, c/o Saville Electronics, 2613 Joppa Terrace, Baltimore, MD 21234.

**Troubleshooting**

**Problem:** Neither the light nor the clock works.

**Solution 1.** Make sure there is power at the outlet by plugging in a different appliance, or by plugging the clock into another outlet.

**Solution 2.** Unplug the clock from the wall, count to 10 and replug it in. Sometimes the clock's computer chip picks up static charges of electricity and unplugging it removes the static.

**Problem:** The clock is not doing what I expect it to.

**Solution:** Unplug the clock from the wall, count to 10 then replug it in and switch the power on to restore the default settings. Reprogram options you use.

..... **EZ Wake Warranty Form** .....

**NAME:** \_\_\_\_\_

**ADDRESS:** \_\_\_\_\_

**PHONE:** \_\_\_\_\_

**DATE PURCHASED:** \_\_\_\_\_

**PLACE PURCHASED:** \_\_\_\_\_

**MODEL 88620-SERIAL#:** \_\_\_\_\_

# *Instructions for the* **EZ WAKE** <sup>TM</sup>

## *SunRise Clock*

### **Bio-Brite Inc.**

4340 East-West Hwy. Suite 401S

Bethesda, MD 20814

301-961-5940

Email: [biobrite@aol.com](mailto:biobrite@aol.com)

Website [www.biobrite.com](http://www.biobrite.com)

